

## 兔抗 BPIFB1 多克隆抗体

中文名称：兔抗 BPIFB1 多克隆抗体

英文名称： Anti-BPIFB1 rabbit polyclonal antibody

别名： BPI fold containing family B member 1; LPLUNC1; C20orf114

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： BPIFB1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BPIFB1
<b>Immunogen:</b>	Fusion protein of human BPIFB1
<b>Full name:</b>	BPI fold containing family B member 1
<b>Synonyms:</b>	LPLUNC1; C20orf114

<b>SwissProt:</b>	Q8TDL5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	30-150

