

## 兔抗 BPIFB2 多克隆抗体

中文名称：兔抗 BPIFB2 多克隆抗体

英文名称： Anti-BPIFB2 rabbit polyclonal antibody

别名： RYSR, BPIL1, LPLUNC2, C20orf184, dJ726C3.2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： BPIFB2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the lipid transfer/lipo polysaccharide binding protein (LT/LBP) gene family. It is highly expressed in hypertrophic tonsils. This gene and three other members of the LT/LBP gene family form a cluster on the long arm of chromosome 20.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BPIFB2
<b>Immunogen:</b>	Synthetic peptide of human BPIFB2
<b>Full name:</b>	BPI fold containing family B, member 2
<b>Synonyms :</b>	RYSR, BPIL1, LPLUNC2, C20orf184, dJ726C3.2
<b>SwissProt:</b>	Q8N4F0

<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human liver cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

