

兔抗 BPI 多克隆抗体

中文名称：兔抗 BPI 多克隆抗体

英文名称： Anti-BPI rabbit polyclonal antibody

别名： rBPI; BPIFD1

相关类别： 一抗

抗原： BPI

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a lipopolysaccharide binding protein. It is associated with human neutrophil granules and has bactericidal activity on gram-negative organisms. The cytotoxic action of BPI is limited to many species of Gram-negative bacteria; this specificity may be explained by a strong affinity of the very basic N-terminal half for the negatively charged lipopolysaccharides that are unique to the Gram-negative bacterial outer envelope. Has antibacterial activity against the Gram-negative bacterium *P.aeruginosa*, this activity is inhibited by LPS from *P.aeruginosa*.

Applications:	ELISA, IHC
Name of antibody:	BPI
Immunogen:	Fusion protein of human BPI
Full name:	bactericidal/permeability-increasing protein
Synonyms :	rBPI; BPIFD1
SwissProt:	P17213
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	50-200

