

RAB11A 抗原（重组蛋白）

中文名称：RAB11A 抗原（重组蛋白）

英文名称：RAB11A Antigen (Recombinant Protein)

别名：YL8

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	RAB11A, member RAS oncogene family
Synonyms:	YL8
Swissprot:	P62491
Gene Accession:	BC013348
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene.