

## 兔抗 CYTIP 多克隆抗体

中文名称：兔抗 CYTIP 多克隆抗体

英文名称： Anti-CYTIP rabbit polyclonal antibody

别名： cytohesin 1 interacting protein; HE; B3-1; CASP; CYBR; CYTHIP; PSCDBP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CYTIP

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene contains 2 leucine zipper domains and a putative C-terminal nuclear targeting signal, but does not have any hydrophobic regions. This protein is expressed weakly in resting NK and T cells. The encoded protein modulates the activation of ARF genes by CYTH1. This protein interacts with CYTH1 and SNX27 proteins and may act to sequester CYTH1 protein in the cytoplasm.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CYTIP
<b>Immunogen:</b>	Fusion protein of human CYTIP

<b>Full name:</b>	cytohesin 1 interacting protein
<b>Synonyms:</b>	HE; B3-1; CASP; CYBR; CYTHIP; PSCDBP
<b>SwissProt:</b>	O60759
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human brain
<b>IHC Recommend dilution:</b>	50-100



