

兔抗 CYR1 多克隆抗体

中文名称：兔抗 CYR1 多克隆抗体

英文名称：Anti-CYR1 rabbit polyclonal antibody

别名：cysteine and tyrosine rich 1; C21orf95

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CYR1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein is unknown.
Applications:	ELISA, IHC
Name of antibody:	CYR1
Immunogen:	Fusion protein of human CYR1
Full name:	cysteine and tyrosine rich 1
Synonyms:	C21orf95
SwissProt:	Q96J86
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-300

