

RAB39B 抗原（重组蛋白）

中文名称：RAB39B 抗原（重组蛋白）

英文名称：RAB39B Antigen (Recombinant Protein)

别名：WSMN; MRX72

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	RAB39B, member RAS oncogene family
Synonyms:	WSMN; MRX72
Swissprot:	Q96DA2
Gene Accession:	BC009714
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the Rab family of proteins. Rab proteins are small GTPases that are involved in vesicular trafficking. Mutations in this gene are associated with X-linked mental retardation.