

兔抗 BRAF(Phospho-Ser446) 多克隆抗体

中文名称：兔抗 BRAF(Phospho-Ser446) 多克隆抗体

英文名称： Anti-BRAF(Phospho-Ser446) rabbit polyclonal antibody

别名： NS7; BRAF1; RAFB1; B-RAF1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： BRAF(Phospho-Ser446)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein belonging to the raf/mil family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene are associated with cardiofaciocutaneous syndrome, a disease characterized by heart defects, mental retardation and a distinctive facial appearance.
Applications:	WB
Name of antibody:	BRAF(Phospho-Ser446)
Immunogen:	Peptide sequence around phosphorylation site of Serine 446 (R-D-S(p)-S-D) derived from Human B-RAF.

Full name:	B-Raf proto-oncogene, serine/threonine kinase
Synonyms :	NS7; BRAF1; RAFB1; B-RAF1
SwissProt:	P15056
WB Predicted band size:	84 kDa
WB Positive control:	293 cells lysates treated with EGF
WB Recommended dilution:	500-1000

