

兔抗 DAB2IP 多克隆抗体

中文名称：兔抗 DAB2IP 多克隆抗体

英文名称：Anti-DAB2IP rabbit polyclonal antibody

别名：AIP1; AIP-1; AF9Q34; DIP1/2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DAB2IP

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | DAB2IP is a Ras (MIM 190020) GTPase-activating protein (GAP) that acts as a tumor suppressor. The DAB2IP gene is inactivated by methylation in prostate and breast cancers |
| Applications: | ELISA, IHC |
| Name of antibody: | DAB2IP |
| Immunogen: | Synthetic peptide of human DAB2IP |
| Full name: | DAB2 interacting protein |
| Synonyms : | AIP1; AIP-1; AF9Q34; DIP1/2 |
| SwissProt: | Q5VWQ8 |
| ELISA Recommended dilution: | 1000-2000 |

| | |
|--------------------------------|---|
| IHC positive control: | Human thyroid cancer and Human gastric cancer |
| IHC Recommend dilution: | 25-100 |

