

RAB4A 抗原（重组蛋白）

中文名称：RAB4A 抗原（重组蛋白）

英文名称：RAB4A Antigen (Recombinant Protein)

别名：RAB4; HRES1; HRES-1; HRES-1/RAB4

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	RAB4A, member RAS oncogene family
Synonyms:	RAB4; HRES1; HRES-1; HRES-1/RAB4
Swissprot:	P20338
Gene Accession:	BC002438
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.