

RERG 抗原（重组蛋白）

中文名称：RERG 抗原（重组蛋白）

英文名称：RERG Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

| | |
|---------------------------|--|
| Full name: | RAS-like, estrogen-regulated, growth inhibitor |
| Swissprot: | Q96A58 |
| Gene Accession: | BC007997 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | RERG, a member of the RAS superfamily of GTPases, inhibits cell proliferation and tumor formation (Finlin et al., 2001 [PubMed 11533059]). |