

## 兔抗 DAK 多克隆抗体

中文名称：兔抗 DAK 多克隆抗体

英文名称：Anti-DAK rabbit polyclonal antibody

别名：NET45

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：DAK

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene is a member of the family of dihydroxyacetone kinases, which have a protein structure distinct from other kinases. The product of this gene phosphorylates dihydroxyacetone, and also catalyzes the formation of riboflavin 4',5'-phosphate (aka cyclin FMN) from FAD. Several alternatively spliced transcript variants have been identified, but the full-length nature of only one has been determined.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	DAK
<b>Immunogen:</b>	Synthesized peptide derived from internal of human

	DAK.
<b>Full name:</b>	dihydroxyacetone kinase 2 homolog
<b>Synonyms :</b>	NET45
<b>SwissProt:</b>	Q3LXA3
<b>WB Predicted band size:</b>	59 kDa
<b>WB Positive control:</b>	K562 cells and A549 cells lysates
<b>WB Recommended dilution:</b>	500-3000

