

## RILPL2 抗原（重组蛋白）

中文名称：RILPL2 抗原（重组蛋白）

英文名称：RILPL2 Antigen (Recombinant Protein)

别名：RLP2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	Rab interacting lysosomal protein-like 2
<b>Synonyms:</b>	RLP2
<b>Swissprot:</b>	Q969X0
<b>Gene Accession:</b>	BC013042
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that contains a rab-interacting lysosomal protein-like domain. This protein may be involved in regulating lysosome morphology. This protein may also be a target for the Hepatitis C virus and assist in viral replication. Alternate splicing results in multiple transcript variants.