

## 兔抗 Dap3 多克隆抗体

中文名称：兔抗 Dap3 多克隆抗体

英文名称： Anti-Dap3 rabbit polyclonal antibody

别名： DAP-3, 4921514D13Rik

抗原： Dap3

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

28S ribosomal protein S29, mitochondrial is a protein that in humans is encoded by the DAP3 gene. Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. This protein potentially binds ATP/GTP and might be a functional partner of the mitoribosomal protein S27. Splice variants that differ in the 5' UTR have been found for this gene; both variants encode the same protein. Pseudogenes corresponding to this gene are found on ch

	romosomes 1q and 2q. DAP3 has been shown to interact with FADD, TNFRSF10A, Glucocorticoid receptor and Heat shock protein 90kDa alpha (cytosolic), member A1.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	Dap3
<b>Immunogen:</b>	Synthetic peptide of mouse Dap3
<b>Full name:</b>	death associated protein 3
<b>Synonyms :</b>	DAP-3, 4921514D13Rik
<b>SwissProt:</b>	Q9ER88
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	10-50
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	Hepg2 and Jurkat cells
<b>WB Recommended dilution:</b>	500-1000



