

兔抗 DCAF12 多克隆抗体

中文名称：兔抗 DCAF12 多克隆抗体

英文名称： Anti-DCAF12 rabbit polyclonal antibody

别名： DDB1 and CUL4 associated factor 12; CT102; TCC52; WDR40A; KIAA1892

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DCAF12

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates.
Applications:	ELISA, IHC
Name of antibody:	DCAF12
Immunogen:	Synthetic peptide of human DCAF12
Full name:	DDB1 and CUL4 associated factor 12
Synonyms:	CT102; TCC52; WDR40A; KIAA1892
SwissProt:	Q5T6F0
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer and Human tonsil
IHC Recommend dilution:	50-200



