

RSU1 抗原（重组蛋白）

中文名称：RSU1 抗原（重组蛋白）

英文名称：RSU1 Antigen (Recombinant Protein)

别名：Ras suppressor protein 1; RSP-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 2-277AA amino acids of human RSU1

技术规格

Full name:	Ras suppressor protein 1
Synonyms:	RSP-1
Swissprot:	Q15404
Gene Accession:	BC005993
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that is involved in the Ras signal transduction pathway, growth inhibition, and nerve-growth factor induced differentiation processes, as determined in mouse and human cell line studies. In mouse, the encoded protein was initially isolated based on its ability to inhibit v-Ras transformation. Multiple alternatively spliced transcript variants for this gene have been reported; one of these variants was found only in glioma tumors.