

## 兔抗 BTK(Phospho-Tyr551)多克隆抗体

中文名称：兔抗 BTK(Phospho-Tyr551)多克隆抗体

英文名称：Anti-BTK(Phospho-Tyr551) rabbit polyclonal antibody

别名：Bruton tyrosine kinase; AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1

相关类别：一抗

抗原：BTK(Phospho-Tyr551)

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene plays a crucial role in B-cell development. Mutations in this gene cause X-linked agammaglobulinemia type 1, which is an immunodeficiency characterized by the failure to produce mature B lymphocytes, and associated with a failure of Ig heavy chain rearrangement. Alternative splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	BTK(Phospho-Tyr551)
<b>Immunogen:</b>	Synthetic peptide of human BTK

<b>Full name:</b>	Bruton tyrosine kinase
<b>Synonyms:</b>	AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1
<b>SwissProt:</b>	Q06187
<b>WB Predicted band size:</b>	76 kDa
<b>WB Positive control:</b>	Hela cells
<b>WB Recommended dilution:</b>	500-1000

