

兔抗 DCAF13 多克隆抗体

中文名称: 兔抗 DCAF13 多克隆抗体

英文名称: Anti-DCAF13 rabbit polyclonal antibody

别名: GM83; WDSOF1; HSPC064

相关类别: 一抗

抗原: DCAF13

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	DCAF13 localizes to the granular component of the nucleolus, a region associated with late ribosomal RNA-processing events and assembly of ribosomal particles. Possible role in ribosomal RNA processing (By similarity). May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.
Applications:	ELISA, WB
Name of antibody:	DCAF13
Immunogen:	Synthetic peptide of human DCAF13
Full name:	DDB1 and CUL4 associated factor 13
Synonyms:	GM83; WDSOF1; HSPC064

SwissProt:	Q9NV06
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	51 kDa
WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	200-1000

