

兔抗 BTK 多克隆抗体

中文名称：兔抗 BTK 多克隆抗体

英文名称：Anti-BTK rabbit polyclonal antibody

别名：AT; ATK; BPK; XLA; IMD1; AGMX1; PSCTK1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BTK

标记物：Unconjugate 克

隆类型：rabbit polyclonal

反应种属：Human, Mouse

技术规格

Background:	The protein encoded by this gene plays a crucial role in B-cell development. Mutations in this gene cause X-linked agammaglobulinemia type 1, which is an immunodeficiency characterized by the failure to produce mature B lymphocytes, and associated with a failure of Ig heavy chain rearrangement. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	BTK
Immunogen:	Fusion protein of human BTK
Full name:	Bruton tyrosine kinase

Synonyms:	AT; ATK; BPK; XLA; IMD1; AGMX1; PSCTK1
SwissProt:	Q06187
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human brain
IHC Recommend dilution:	25-100



