

兔抗 BTLA 多克隆抗体

中文名称：兔抗 BTLA 多克隆抗体

英文名称： Anti-BTLA rabbit polyclonal antibody

别名： BTLA1; CD272

相关类别： 一抗

抗原： BTLA

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.
Applications:	ELISA, WB, IHC
Name of antibody:	BTLA
Immunogen:	Fusion protein of human BTLA

Full name:	B and T lymphocyte associated
Synonyms:	BTLA1; CD272
SwissProt:	Q7Z6A9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human gastric cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	33 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	200-1000



