

## 兔抗 DCAF17 多克隆抗体

中文名称：兔抗 DCAF17 多克隆抗体

英文名称：Anti-DCAF17 rabbit polyclonal antibody

别名：DDB1 and CUL4 associated factor 17 C2orf37

相关类别：一抗

抗原：DCAF17

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a nuclear transmembrane protein that associates with cullin 4A/damaged DNA binding protein 1 ubiquitin ligase complex. Mutations in this gene are associated with Woodhouse-Sakati syndrome. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DCAF17
<b>Immunogen:</b>	Synthetic peptide of human DCAF17
<b>Full name:</b>	DDB1 and CUL4 associated factor 17
<b>Synonyms:</b>	C2orf37
<b>SwissProt:</b>	Q5H9S7

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	20-100



