

## S100A14 抗原（重组蛋白）

中文名称：S100A14 抗原（重组蛋白）

英文名称：S100A14 Antigen (Recombinant Protein)

别名：BCMP84; S100A15

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	S100 calcium binding protein A14
<b>Synonyms:</b>	BCMP84; S100A15
<b>Swissprot:</b>	Q9HCY8
<b>Gene Accession:</b>	BC005019
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the S100 protein family which contains an EF-hand motif and binds calcium. The gene is located in a cluster of S100 genes on chromosome 1. Levels of the encoded protein have been found to be lower in cancerous tissue and associated with metastasis suggesting a tumor suppressor function