

## 兔抗 DCBLD2 多克隆抗体

中文名称：兔抗 DCBLD2 多克隆抗体

英文名称： Anti-DCBLD2 rabbit polyclonal antibody

别名： ESDN; CLCP1

相关类别： 一抗

抗原： DCBLD2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

CLCP1 (CUB, LCCL and coagulation factor V/VIII-homology domains protein 1), also known as DCBLD2 (discoïdin, CUB and LCCL domain containing 2) or ESDN, is a 775 amino acid single-pass type I membrane protein that contains one CUB domain, one LCCL domain and one F5/8 type C domain. Expressed at high levels in heart, testis and skeletal muscle, CLCP1 is thought to regulate vascular smooth muscle cell (VSMC) proliferation and remodeling and may be involved in the transformation and metastasis of various cancers, such as metastatic lung cancer and gastric car

	cinoma.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DCBLD2
<b>Immunogen:</b>	Fusion protein of human DCBLD2
<b>Full name:</b>	discoidin, CUB and LCCL domain containing 2
<b>Synonyms :</b>	ESDN; CLCP1
<b>SwissProt:</b>	Q96PD2
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	100-300

