

兔抗 BUB3(Phospho-Tyr207) 多克隆抗体

中文名称：兔抗 BUB3(Phospho-Tyr207) 多克隆抗体

英文名称：Anti-BUB3(Phospho-Tyr207) rabbit polyclonal antibody

别名：BUB3L; hBUB3

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：BUB3(Phospho-Tyr207)

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Applications:	WB
Name of antibody:	BUB3(Phospho-Tyr207)
Immunogen:	Peptide sequence around phosphorylation site of Tyrosine 207 (V-E-Y(p)-L-D) derived from Human Bub3.
Full name:	BUB3 mitotic checkpoint protein
Synonyms :	BUB3L; hBUB3

SwissProt:	O43684
WB Predicted band size:	37 kDa
WB Positive control:	U87 and U87 EGFRvIII cells in Mitosis lysates
WB Recommended dilution:	500-1000

