

兔抗 C1QBP 多克隆抗体

中文名称：兔抗 C1QBP 多克隆抗体

英文名称： Anti-C1QBP rabbit polyclonal antibody

别名： p32; HABP1; gC1qR; GC1QBP; SF2p32; gC1Q-R

抗原： C1QBP

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The human complement subcomponent C1q associates with C1r and C1s in order to yield the first component of the serum complement system. The protein encoded by this gene is known to bind to the globular heads of C1q molecules and inhibit C1 activation. This protein has also been identified as the p32 subunit of pre-mRNA splicing factor SF2, as well as a hyaluronic acid-binding protein.
Applications:	WB, IHC
Name of antibody:	C1QBP
Immunogen:	Synthesized peptide derived from internal of human C

	1QBP.
Full name:	complement component 1, q subcomponent binding p rotein
Synonyms :	p32; HABP1; gC1qR; GC1QBP; SF2p32; gC1Q-R
SwissProt:	Q07021
IHC positive control:	Human tonsil tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	31 kDa
WB Positive control:	Jurkat cells, K562 cells and Huvec cells lysates
WB Recommended dilution:	500-3000



