

## 兔抗 C1QL2 多克隆抗体

中文名称：兔抗 C1QL2 多克隆抗体

英文名称： Anti-C1QL2 rabbit polyclonal antibody

别名： CTRP10; C1QTNF10

抗原： C1QL2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	May regulate the number of excitatory synapses that are formed on hippocampus neurons. Has no effect on inhibitory synapses By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	C1QL2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human C1QL2.
<b>Full name:</b>	complement component 1, q subcomponent-like 2
<b>Synonyms :</b>	CTRP10; C1QTNF10
<b>SwissProt:</b>	Q7Z5L3
<b>WB Predicted band size:</b>	29 kDa

<b>WB Positive control:</b>	Huvec cells and Hela cells lysates
<b>WB Recommended dilution:</b>	500-3000

