

## 兔抗 DCTD 多克隆抗体

中文名称：兔抗 DCTD 多克隆抗体

英文名称： Anti-DCTD rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： DCTD

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene catalyzes the deamination of dCMP to dUMP, the nucleotide substrate for thymidylate synthase. The encoded protein is allosterically activated by dCTP and inhibited by dTTP, and is found as a homohexamer. This protein uses zinc as a cofactor for its activity. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DCTD
<b>Immunogen:</b>	Full length fusion protein(178FULL)
<b>Full name:</b>	dCMP deaminase
<b>SwissProt:</b>	P32321
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC Positive control:</b>	Human liver cancer;
<b>IHC Recommended dilution:</b>	25-100

