

兔抗 C1QTNF6 多克隆抗体

中文名称：兔抗 C1QTNF6 多克隆抗体

英文名称： Anti-C1QTNF6 rabbit polyclonal antibody

别名： CTFP6; CTRP6; ZACRP6

抗原： C1QTNF6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein is not uniformly distributed but rather focally localized in hepatocellular carcinoma cells and the tumor endothelial cells. It may play a role in metabolism and immunity. CTRPs acts as new molecular mediators connecting inflammatory and metabolic diseases.
Applications:	ELISA, IHC
Name of antibody:	C1QTNF6
Immunogen:	Fusion protein of human C1QTNF6
Full name:	C1q and tumor necrosis factor related protein 6
Synonyms :	CTFP6; CTRP6; ZACRP6
SwissProt:	Q9BXI9

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100

