

SERTAD1 抗原（重组蛋白）

中文名称：SERTAD1 抗原（重组蛋白）

英文名称： SERTAD1 Antigen (Recombinant Protein)

别名： SEI1; TRIP-Br1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	SERTA domain containing 1
Synonyms:	SEI1; TRIP-Br1
Swissprot:	Q9UHV2
Gene Accession:	BC002670
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Acts at E2F-responsive promoters as coregulator to integrate signals provided by PHD- and/or bromodomain-containing transcription factors. Stimulates E2F1/TFDP1 transcriptional activity. Renders the activity of cyclin D1/CDK4 resistant to the inhibitory effects of CDKN2A/p16INK4A.