

## SFN 抗原（重组蛋白）

中文名称：SFN 抗原（重组蛋白）

英文名称： SFN Antigen (Recombinant Protein)

别名： YWHAS

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	stratifin
<b>Synonyms:</b>	YWHAS
<b>Swissprot:</b>	P31947
<b>Gene Accession:</b>	BC000329
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. When bound to KRT17, regulates protein synthesis and epithelial cell growth by stimulating Akt/mTOR pathway. May also regulate MDM2 autoubiquitination and degradation and ther

eby activate p53/TP53.