

## 兔抗 DCTN5 多克隆抗体

中文名称：兔抗 DCTN5 多克隆抗体

英文名称： Anti-DCTN5 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： DCTN5

宿 主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a subunit of dynactin, a component of the cytoplasmic dynein motor machinery involved in minus-end-directed transport. The encoded protein is a component of the pointed-end subcomplex and is thought to bind membranous cargo. A pseudogene of this gene is located on the long arm of chromosome 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DCTN5
<b>Immunogen:</b>	Fusion protein of human DCTN5
<b>Full name:</b>	dynactin 5 (p25)
<b>SwissProt:</b>	Q9BTE1

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	25-100

