

兔抗 DCTN6 多克隆抗体

中文名称：兔抗 DCTN6 多克隆抗体

英文名称： Anti-DCTN6 rabbit polyclonal antibody

别名： WS3; p27; WS-3

相关类别： 一抗

抗原： DCTN6

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains an RGD (Arg-Gly-Asp) motif in the N-terminal region, which confers adhesive properties to macromolecular proteins like fibronectin. It shares a high degree of sequence similarity with the mouse homolog, which has been suggested to play a role in mitochondrial biogenesis. The exact biological function of this gene is not known.
Applications:	ELISA, WB, IHC
Name of antibody:	DCTN6
Immunogen:	Fusion protein of human DCTN6
Full name:	dynactin 6

Synonyms:	WS3; p27; WS-3
SwissProt:	O00399
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	21 kDa
WB Positive control:	A549 cells
WB Recommended dilution:	200-1000



