

兔抗 C1R 多克隆抗体

中文名称： 兔抗 C1R 多克隆抗体

英文名称： Anti-C1R rabbit polyclonal antibody

别名： complement C1r subcomponent;

相关类别： 一抗

抗原： C1R

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	C1r B chain is a serine protease that combines with C1q and C1s to form C1, the first component of the classical pathway of the complement system.
Applications:	ELISA, WB, IHC
Name of antibody:	C1R
Immunogen:	Fusion protein of human C1R
Full name:	complement C1r subcomponent
SwissProt:	P00736
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;
IHC Recommend dilution:	40-200

WB Predicted band size:	80 kDa
WB Positive control:	Human plasma solution
WB Recommended dilution:	500-2000



