

兔抗 C2CD2 多克隆抗体

| | |
|-------|--|
| 中文名称: | 兔抗 C2CD2 多克隆抗体 |
| 英文名称: | Anti-C2CD2 rabbit polyclonal antibody |
| 别名: | C2 calcium-dependent domain containing 2; TMEM24L; C21orf25; C21orf258 |
| 相关类别: | 一抗 |
| 储存: | 冷冻 (-20℃) |
| 宿主: | Rabbit |
| 抗原: | C2CD2 |
| 反应种属: | Human |
| 标记物: | Unconjugate |
| 克隆类型: | rabbit polyclonal |

技术规格

| | |
|------------------------------------|--|
| Background: | C2CD2 is an 696 amino acid protien that is secreted. The C2CD2 gene product has been provisionally designated C 2CD2 pending further characterization. |
| Applications: | ELISA, IHC |
| Name of antibody: | C2CD2 |
| Immunogen: | Fusion protein of human C2CD2 |
| Full name: | C2 calcium-dependent domain containing 2 |
| Synonyms: | TMEM24L; C21orf25; C21orf258 |
| SwissProt: | Q9Y426 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human liver cancer |

IHC Recommend dilution: 25-100

