

## 兔抗 C2CD6 多克隆抗体

中文名称：兔抗 C2CD6 多克隆抗体

英文名称：Anti-C2CD6 rabbit polyclonal antibody

别名：C2 calcium dependent domain containing 6; ALS2CR11

相关类别：一抗

抗原：C2CD6

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	An autosomal recessive form of juvenile amyotrophic lateral sclerosis was originally mapped to a region of chromosome 2 that includes this gene. The encoded protein contains a calcium-dependent membrane targeting C2 domain. This domain is often found in proteins that are involved in membrane trafficking and signal transduction.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	C2CD6
<b>Immunogen:</b>	Fusion protein of human C2CD6
<b>Full name:</b>	C2 calcium dependent domain containing 6

<b>Synonyms:</b>	ALS2CR11
<b>SwissProt:</b>	Q53TS8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

