

SMS 抗原（重组蛋白）

中文名称：SMS 抗原（重组蛋白）

英文名称：SMS Antigen (Recombinant Protein)

别名：spermine synthase; SRS; SpS; MRSR; SPMSY

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 167-366 amino acids of human SMS

技术规格

Full name:	spermine synthase
Synonyms:	SRS; SpS; MRSR; SPMSY
Swissprot:	P52788
Gene Accession:	BC009898
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein belonging to the spermidine/spermin synthase family. Pseudogenes of this gene are located on chromosomes 1, 5, 6 and X. Mutations in this gene are associated with X-linked Snyder-Robinson mental retardation syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.