

兔抗 C3orf38 多克隆抗体

中文名称：兔抗 C3orf38 多克隆抗体

英文名称：Anti-C3orf38 rabbit polyclonal antibody

别名：chromosome 3 open reading frame 38

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：C3orf38

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in apoptosis regulation.
Applications:	ELISA, IHC
Name of antibody:	C3orf38
Immunogen:	Fusion protein of human C3orf38
Full name:	chromosome 3 open reading frame 38
SwissProt:	Q5JPI3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	100-300

