

兔抗 C4orf46 多克隆抗体

中文名称：兔抗 C4orf46 多克隆抗体

英文名称：Anti-C4orf46 rabbit polyclonal antibody

别名：RCDG1

抗原：C4orf46

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a small, conserved protein of unknown function that is expressed in a variety of tissues. There are pseudogenes for this gene on chromosomes 6, 8, 16, and X. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	C4orf46
Immunogen:	Full length fusion protein
Full name:	chromosome 4 open reading frame 46
Synonyms :	RCDG1
SwissProt:	Q504U0

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

