

SNX11 抗原（重组蛋白）

中文名称：SNX11 抗原（重组蛋白）

英文名称：SNX11 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	sorting nexin 11
Swissprot:	Q9Y5W9
Gene Accession:	BC000768
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members. This gene encodes a protein of unknown function. This gene results in two transcript variants differing in the 5' UTR, but encoding the same protein.