

SSX1 抗原（重组蛋白）

中文名称：SSX1 抗原（重组蛋白）

英文名称：SSX1 Antigen (Recombinant Protein)

别名：SSRC; CT5.1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	synovial sarcoma, X breakpoint 1
Synonyms:	SSRC; CT5.1
Swissprot:	Q16384
Gene Accession:	BC128611
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy.