

兔抗 C9 多克隆抗体

中文名称：兔抗 C9 多克隆抗体

英文名称： Anti-C9 rabbit polyclonal antibody

别名： C9D; ARMD15

储存： 冷冻（-20℃） 避光

抗原： C9

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes the final component of the complement system. It participates in the formation of the Membrane Attack Complex (MAC). The MAC assembles on bacterial membranes to form a pore, permitting disruption of bacterial membrane organization. Mutations in this gene cause component C9 deficiency.
Applications:	WB
Name of antibody:	C9
Immunogen:	Synthesized peptide derived from internal of human C9.
Full name:	complement component 9

Synonyms :	C9D; ARMD15
SwissProt:	P02748
WB Predicted band size:	63 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000

