

兔抗 DDIT3(Phospho-Ser30) 多克隆抗体

中文名称：兔抗 DDIT3(Phospho-Ser30) 多克隆抗体

英文名称： Anti-DDIT3(Phospho-Ser30) rabbit polyclonal antibody

别名： CHOP; CEBPZ; CHOP10; CHOP-10; GADD153

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： DDIT3(Phospho-Ser30)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the CCAAT/enhancer-binding protein (C/EBP) family of transcription factors. The protein functions as a dominant-negative inhibitor by forming heterodimers with other C/EBP members, such as C/EBP and LAP (liver activator protein), and preventing their DNA binding activity. The protein is implicated in adipogenesis and erythropoiesis, is activated by endoplasmic reticulum stress, and promotes apoptosis.
Applications:	WB
Name of antibody:	DDIT3(Phospho-Ser30)
Immunogen:	Peptide sequence around phosphorylation site of serine30

	V-L-S(p)-S-D) derived from Human CHOP .
Full name:	DNA-damage-inducible transcript 3
Synonyms :	CHOP; CEBPZ; CHOP10; CHOP-10; GADD153
SwissProt:	P35638
WB Predicted band size:	19 kDa
WB Positive control:	Jurkat cells lysates treated with PMA
WB Recommended dilution:	500-1000

