

## STX11 抗原（重组蛋白）

中文名称：STX11 抗原（重组蛋白）

英文名称：STX11 Antigen (Recombinant Protein)

别名：FHL4; HLH4; HPLH4

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	syntaxin 11
<b>Synonyms:</b>	FHL4; HLH4; HPLH4
<b>Swissprot:</b>	O75558
<b>Gene Accession:</b>	BC033519
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the syntaxin family. Syntaxins have been implicated in the targeting and fusion of intracellular transport vesicles. This family member may regulate protein transport among late endosomes and the trans-Golgi network. Mutations in this gene have been associated with familial hemophagocytic lymphohistiocytosis.