

## 兔抗 C14orf180 多克隆抗体

- 中文名称：兔抗 C14orf180 多克隆抗体
- 英文名称：Anti-C14orf180 rabbit polyclonal antibody
- 别名：chromosome 14 open reading frame 180; NRAC; C14orf77
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：C14orf180
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	C14orf180 is a multi-pass membrane protein. The function of the C14orf180 protein is not known.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	C14orf180
<b>Immunogen:</b>	Synthetic peptide of human C14orf180
<b>Full name:</b>	chromosome 14 open reading frame 180
<b>Synonyms:</b>	NRAC; C14orf77
<b>SwissProt:</b>	Q8N912
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100

