

兔抗 C19orf57 多克隆抗体

中文名称：兔抗 C19orf57 多克隆抗体

英文名称：Anti-C19orf57 rabbit polyclonal antibody

别名：chromosome 19 open reading frame 57

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：C19orf57

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, WB, IHC
Name of antibody:	C19orf57
Immunogen:	Fusion protein of human C19orf57
Full name:	chromosome 19 open reading frame 57
SwissProt:	Q0VDD7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	70 kDa
WB Positive control:	PC3 cell lysate

WB Recommended dilution: 200-1000

