

TAC1 抗原（重组蛋白）

中文名称：TAC1 抗原（重组蛋白）

英文名称：TAC1 Antigen (Recombinant Protein)

别名：NK2; NPK; NKNA; TAC2; Hs.2563

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	tachykinin, precursor 1
Synonyms:	NK2; NPK; NKNA; TAC2; Hs.2563
Swissprot:	P20366
Gene Accession:	BC018047
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues. Multiple transcript variants encoding different isoforms have been f

ound for this gene.