

TAGLN 抗原（重组蛋白）

中文名称：TAGLN 抗原（重组蛋白）

英文名称：TAGLN Antigen (Recombinant Protein)

别名：SM22; SMCC; TAGLN1; WS3-10

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	transgelin
Synonyms:	SM22; SMCC; TAGLN1; WS3-10
Swissprot:	Q01995
Gene Accession:	BC004927
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a transformation and shape-change sensitive actin cross-linking/gelling protein found in fibroblasts and smooth muscle. Its expression is down-regulated in many cell lines, and this down-regulation may be an early and sensitive marker for the onset of transformation. A functional role of this protein is unclear. Two transcript variants encoding the same protein have been found for this gene.