

TCEAL4 抗原（重组蛋白）

中文名称：TCEAL4 抗原（重组蛋白）

英文名称：TCEAL4 Antigen (Recombinant Protein)

别名：WEX7; NPD017

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	transcription elongation factor A (SII)-like 4
Synonyms:	WEX7; NPD017
Swissprot:	Q96E15
Gene Accession:	BC012296
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. This family is comprised of nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X c chromosome. Alternatively splicing results in multiple transcript variants. There is a pseudogene for this gene on chromosome 13.