

TFDP2 抗原（重组蛋白）

中文名称：TFDP2 抗原（重组蛋白）

英文名称：TFDP2 Antigen (Recombinant Protein)

别名：DP2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	transcription factor Dp-2 (E2F dimerization partner 2)
Synonyms:	DP2
Swissprot:	Q14188
Gene Accession:	BC021113
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The gene is a member of the transcription factor DP family. The encoded protein forms heterodimers with the E2F transcription factors resulting in transcriptional activation of cell cycle regulated genes. Alternative splicing results in multiple transcript variants.