

## TGIF2LX 抗原（重组蛋白）

中文名称：TGIF2LX 抗原（重组蛋白）

英文名称：TGIF2LX Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 236 amino acids of human TGIF2LX

技术规格

<b>Full name:</b>	TGFB-induced factor homeobox 2-like, X-linked
<b>Synonyms:</b>	TGIFLX
<b>Swissprot:</b>	Q8IUE1
<b>Gene Accession:</b>	BC117398
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. Testis-specific expression suggests that this gene may play a role in spermatogenesis. A homolog of this gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition.